

OEM Pump Series

Four-Piston Diaphragm Pumps



QUATTROFLOW® OEM PUMP SERIES FOR IMPROVED SYSTEM INTEGRATION

The OEM Pump Series is an integrable (panel-mount) pump series without housing. It comes with small footprint and mounting solutions designed for improved system integration. The attached synchronous servo motor allows a wide flow range and multiple connections to external PLC systems.

OEM SERIES


- For Quattroflow® multiple- and single-use pumps
- Synchronous servomotor option in 230V / 400V for extended flow ranges
- SEW servo motor as standard
- OEM Pump flange option prepared for multiple AllenBradley motor types
- Available for Quattroflow pump models QF30, QF150, QF1200, QF2500, QF4400 and QF5k

ACCESSORIES:


- VFD with optional communication cards (ProfiNET, EtherNET IP, ...) for connection to external PLC system
- Applicable cables in multiple length options
- Gasket solutions to seal pump and system connection point

Pump Size	Panel Mount Sealing	Servo Motor	Motor Cable	Drive Inverter	Communication to PLC
QF30	O-Ring	SEW Allen Bradley*	2m 5m 10m	SEW Allen Bradley*	Analog (4-20mA) ProfiNet EtherNET IP
QF150					
QF1200	Gasket Tunnel				
QF2500					
QF4400					
QF5k					

* Adaptions to several Allen Bradley servo types available.



Quattroflow® OEM Pump Range



* Actual flow-range depends on motor and drive selection



DESIGN FOR LIFE



WARNING: The customer is solely responsible for determining the necessary or most suitable product or service for any particular use case or application. All such risks are assumed by the customer. Specifications subject to change without notice.
WARRANTY: For product warranty information, please contact: info@cpc-bio.com or +1 800-519-7633.
 For trademark information, please visit: www.cpc-bio.com/trademarks-and-patents
 © 2026 CPC BIOTECH. ALL RIGHTS RESERVED.
 CPCBio26001

CPC BIOTECH
 2820 Cleveland Ave. N.
 Roseville, MN 55113
 Phone: +1 800-519-7633
 Email: info@cpc-bio.com

cpc-bio.com